









## Immunofluorescence Labeling Protocol for 2D Cell Culture

STEP		TEMPERATURE	TIME*
1	Warm-up Fixative (4% Paraformaldehyde), PBS**-Cello™-IF3-Primary Antibody Solution-Secondary Antibody Solution ***-Mounting Medium****	37 °C	10 min
2	 Aspirate medium, Fix the sample	37 °C	10 min
3	 Incubate with Cello™-IF2	37 °C	3X, 5 min
4	 Incubate with Primary Antibody in Cello™-IF2	37 °C	30 min-1 hr
5	 Incubate with Cello™-IF	37 °C	3X, 5 min
6	 Incubate with Secondary Antibody in Cello™-IF3, Keep in dark	37 °C	30 min-1 hr
7	 Incubate with PBS, Keep in dark	37 °C	3X, 5 min
8	 Cover with Mounting Medium, Keep in dark	RT	
	 Close the lid tightly, Keep in dark until examination.	4 °C	

**TOTAL : 2.5-4 hr**

\*Timings can be adjusted according to the sample and the antibodies.

\*\*PBS: Phosphate Buffered Saline

\*\*\*Dilute primary and secondary antibodies in Cello-IF.2

\*\*\*\*Mounting Medium:1,2 uL Hoechst or DAPI+675 uL PBS + 675 uL Glycerol

### TIPS:

-Using a cell culture container with high optical quality is highly recommended (e.g., glass-bottomed dishes and well-plates, multi-chambered slides, Cello™ - M). The container here represents Cello™ - M.

-Low-speed orbital shaker can be used to facilitate penetration.

-Most samples are examined on the same day,. Nuclear staining improves after 24 hours in some thick samples.

-Start with the validated dilution rates for each antibody. Adjustment of the dilution rates may be needed based on the nature of the sample and the antibody.

-